

Stand run	Test type	Z156	Z108
1	156 ref	0	9999
2	156 can	0	9999
3	156 can	156	9999
4	156 can	312	9999
5	156 can	468	9999
6	156 can	624	9999
7	156 can	780	9999
8	156 can	936	9999
9	156 can	1092	9999
10	156 can	1248	9999
11	156 can	1404	9999

single-ref, all 156 candidate tests

Stand run	Test type	Z156	Z108
1	156 ref	0	9999
2	108 ref	0	0
3	156 can	0	0
4	156 can	156	156
5	156 can	312	312
6	156 can	468	468
7	156 can	624	624
8	156 can	780	780
9	156 can	936	936
10	156 can	1092	1092
11	156 can	1248	1248
12	156 can	1404	1404

seq dual-ref, all 156 candidate tests

Stand run	Test type	Z156	Z108
1	156 ref	0	9999
2	156 can	0	9999
3	156 can	156	9999
4	156 can	312	9999
5	108 ref	468	0
6	108 can	468	0
7	108 can	576	108
8	108 can	684	216
9	156 can	792	324
10	156 can	948	480
11	156 can	1104	636
12	156 can	1260	792

non-seq dual-ref, 108 & 156 candidate tests